



3M STAINLESS STEEL CROWN/DENOVO BAND

From Dr. LaRee Johnson

Primary 3M SSCs 1st Molar	DENOVO BAND SIZE		Primary 3M SSCs 2nd Molar	DENOVO BAND SIZE		Permanent 3M SSCs	DENOVO BAND SIZE	
	Upper	Lower		Upper	Lower		Upper	Lower
D2	25	25 + crimp	E2	32	31	2	36 1/2	34 1/2
D3	26	25	E3	33 1/2	32	3	37 1/2	36
D4	27	26	E4	34 1/2	33	4	39	37 1/2
D5	28	27	E5	36	34 - 34 1/2	5	40	38 1/2
D6	29	29	E6	37	35 1/2	6	-	39 1/2
D7	30	30	E7	38	36 1/2	7	-	41



NuSmile ZR CROWN/DENOVO BAND

From Dr. LaRee Johnson

Primary NuSmile ZR 1st Molar	DENOVO BAND SIZE				Primary NuSmile ZR 2nd Molar	DENOVO BAND SIZE				Permanent NuSmile ZR	DENOVO BAND SIZE	
	Upper Regular Narrow		Lower Regular Narrow			Upper Regular Narrow		Lower Regular Narrow			Upper	Lower
D1	24 + crimp	23	23 ½	23	E1	28 ½	29	30	28			
D2	25	24	24 ½	24	E2	32	31	32	29	2	37	35
D3	26	25	26	25	E3	34	32	33	30	3	38 ½	36
D4	27	26	27 ½	26 ½	E4	35	33	34	32	4	40	37
D5	28 ½	28	28 ½	27 ½	E5	36	35	35	33	5	41	38 ½
D6	29 ½	29	30	28 ½	E6	38	36+	36	34	6		39 ½
D7	31	30	31 ½	30 ½	E7	39 ½	38	38	35+	7		41

NuSmile Sizes for Band & Loop/Distal Shoe

<i>Maxillary</i>	<i>Mandibular</i>
EUR2 - 31	ELL2 - 30
EUR3 - 32.5	ELL3 - 31
EUR4 - 33.5	ELL4 - 32.5
EUR5 - 35	ELL5 - 33.5
EUR6 - 36	ELL6 - 35
EUR7 - 37.5	ELL7 - 37
DUR2 - 23.5	DLL2 - 24
DUR3 - 24.5	DLL3 - 25.5
DUR4 - 26	DLL4 - 26.5
DUR5 - 28	DLL5 - 27
DUR6 - 28.5	DLL6 - 28.5
DUR7 - 30	DLL7 - 30
DUL2 - 23.5	DLR2 - 24
DUL3 - 24.5	DLR3 - 25.5
DUL4 - 26	DLR4 - 26.5
DUL5 - 28.5	DLR5 - 27
DUL6 - 28.5	DLR6 - 28.5
DUL7 - 30	DLR7 - 30
EUL2 - 31	ELR2 - 30
EUL3 - 32.5	ELR3 - 31
EUL4 - 33.5	ELR4 - 32.5
EUL5 - 35	ELR5 - 33.5
EUL6 - 36	ELR6 - 35
EUL7 - 37.5	ELR7 - 37